

DRY HEAT STERILIZER

Steri-Dent Dry Heat Sterilizers bring high quality, rugged performance, and total dependability to your practice.

The Best Method for Sharp Instruments

Steri-Dent Dry Heat Sterilizers kill bacteria without moisture and will not dull sharp instruments. Plus, they can be used to dry wet packs or as a back-up sterilizer in your practice — the perfect complement to an autoclave.

Safe and Simple to Operate

Dry heat sterilization works without pressure, steam, or chemicals, which increases safety. Steri-Dent products feature a rapid-start heating element for fast operation.

High-Quality

Stainless steel construction means no rust —ever!

Maintenance-Free

Unlike wet heat methods, Steri-Dent Dry Heat Sterilizers require no routine cleaning.

Convenient

Utility trays can be sterilized with lids on and stored without contamination. Steri-Dent Sterilizers also have automatic timers with heat "hold" capability for continuous use in a busy practice.

Value-Priced

Low purchase cost and low operating cost—no distilled water to buy!



Model #200 - (two trays)



Model #300 - (three trays)

STERI-DENT

www.sterident.com

STERI-DENT

Sterilizer Specifications

Model	Tray	Tray Size		Outside dimensions			Inside (useful) space			Shipping Weight	Power Consumption
		W	H	W	H	D	W	H	D		
No. 200	2	10"	6 1/4"	15 1/2"	10 1/2"	9 1/2"	12 1/2"	5"	7 1/2"	17 lbs	500 Watts
No. 300	3	11 3/4"	7 1/4"	19"	14 1/2"	11"	15"	7 1/2"	9"	27 lbs	735 Watts

Operates on 115V. Also available in 220/230V.



Dry Heat Self-Sealing Bags

Use Dry Heat Sterilizing Bags as part of your infection control procedure for an added measure of safety. The azure strip turns brown when it has been exposed to high temperature. The paper used is specifically formulated for dry heat sterilizers. 2.5" x 1.5" x 10.5" (6.4cm x 3.8cm x 26.7cm)

cpac
Imaging